

## Certificate of Quality

**PRODUCT:** PBS(10X), pH7.4

**Cat. No.:** BR110-5L

**Lot No.:**

**COMPONENTS:** 1.37M NaCl  
0.1M Na<sub>2</sub>HPO<sub>4</sub>  
27mM KCl  
18mM KH<sub>2</sub>PO<sub>4</sub>  
pH 7.4

**DESCRIPTION:** Results: Pass

**STORAGE:** Store at room temperature.

**CAUTION:** Avoid contact with eyes, skin and clothing.  
Wash all areas of contact with copious amounts of water.  
For laboratory or further manufacturing use only.