

Certificate of Quality

PRODUCT: RNase A pancreatic
Cat. No.: BR380-500MG
Lot No.:
COMPONENTS: Ribonuclease A
DESCRIPTION: White to light beige powder.

	Min	Max	Units	Disposition
Activity (Kunitz, Protein @37C)	60		u/mg	PASS
Reassay Date	REPORT			PASS
Solubility (1%, Water)	PASS		P/F	PASS

STORAGE: Store frozen. Keep dry. Warm to room temperature before opening.

CAUTION: Avoid contact with eyes, skin and clothing.
Wash all areas of contact with copious amounts of water.
For laboratory or further manufacturing use only.