

Certificate of Quality

PRODUCT: Boric Acid
Cat. No.: BR420-1KG
Lot No.:
COMPONENTS: Boric Acid
DESCRIPTION: ACS GRADE
CAS NO. 10043-35-3
Formula H3BO3
Molecular Weight 61.83

	Min	Max	Units	Disposition
Calcium		0.005	%	PASS
Chloride		0.001	%	PASS
Heavy Metals (as Pb)		0.001	%	PASS
Insolubles (Methanol)		0.005	%	PASS
Iron		0.001	%	PASS
Nonvolatiles (Methanol)		0.05	%	PASS
Phosphate		0.001	%	PASS
Purity	99.5		%	PASS
Sulfate		0.01	%	PASS

STORAGE: Store at room temperature.

CAUTION: Avoid contact with eyes, skin and clothing.
Wash all areas of contact with copious amounts of water.
For laboratory or further manufacturing use only.