

## Certificate of Quality

**PRODUCT:** 10X TBST Buffer

**Cat. No.:** BR510-1L

**Lot No.:**

**COMPONENTS:** 0.5M Tris-Base

1.5M NaCl

0.2% Tween 20

pH 7.6

**DESCRIPTION:** Results: Pass

**STORAGE:** Store at room temperature.

**CAUTION:** Avoid contact with eyes, skin and clothing.  
Wash all areas of contact with copious amounts of water.  
For laboratory or further manufacturing use only.