

## Certificate of Quality

**PRODUCT:** 10X SC Buffer with Tween-20

**Cat. No.:** BR520-1L

**Lot No.:**

**COMPONENTS:** 0.100M Sodium Citrate

0.5% Tween 20

pH 6.0

**DESCRIPTION:** Results: Pass

**STORAGE:** Store at room temperature or at 4°C for longer storage.

**CAUTION:** Avoid contact with eyes, skin and clothing.  
Wash all areas of contact with copious amounts of water.  
For laboratory or further manufacturing use only.