

Western Blot *Ab* Signal Enhancer


Cat. No.: BR570-500

Size: 500 mL


DESCRIPTION: By using Western Blot *Ab* Signal Enhancer as antibody diluents in immuassays, such as western blot and ELISA, sensitivity can be enhanced, background can be reduced, and antigen-antibody reactions can be improved.

Furthermore, Western Blot *Ab* Signal Enhancer is protein-free, so is perfect to work with phosphorylated protein without any interference.


➤ Enhance Sensitivity




➤ Less Antibody Usage



➤ Lower Background/Noises



➤ Reduce Non-Specific Bindings



INSTRUCTIONS:

Simply dilute *primary antibody with Western Blot *Ab* Signal Enhancer to the concentration desired, and then perform assays as usual.

*Based on customers' feedback, secondary antibody diluted by Western Blot *Ab* Signal Enhancer may not work optimally.

STORAGE:

Store at 4°C. Valid for two years from the date of receipt.

For Research Using Only.